

anti- NPHS2 antibody

Product Information

Catalog No.:	FNab05819
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in the regulation of glomerular permeability, acting probably as a linker between the plasma membrane and the cytoskeleton.

Immunogen information

Immunogen:	nephrosis 2, idiopathic, steroid-resistant(podocin)
Synonyms:	NPHS2, PDCN, Podocin, SRN1
Observed MW:	34kd
UniprotID :	Q9NP85

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:200-1:1000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab05819 (NPHS2 antibody) at dilution of 1:500



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab05819(NPHS2 antibody) at dilution of 1:500