

---

## NRF2, NFE2L2 antibody

### Product Information


Catalog No.:	FNab05855
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	nuclear factor(erythroid-derived 2)-like 2
Synonyms:	HEBP1, NF E2 related factor 2, NFE2 related factor 2, NFE2L2, NRF2
Observed MW:	70-100 kDa
Uniprot ID :	Q16236

### Application


Reactivity:	Human, Mouse, Rat, Zebrafish
Tested Application:	ELISA, IP, WB, IHC, FC
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab05855(NRF2, NFE2L2 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



IP Result of anti-NRF2, NFE2L2 (IP: FNab05855, 3ug; Detection: FNab05855 1:300) with mouse kidney tissue lysate 4000ug.



DMSO treated HeLa cells were subjected to SDS PAGE followed by western blot with FNab05855(NRF2, NFE2L2 Antibody) at dilution of 1:600