

## anti- NRSN2 antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab05860   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

May play a role in maintenance and/or transport of vesicles.

### Immunogen information

|              |             |
|--------------|-------------|
| Immunogen:   | neurensin 2 |
| Synonyms:    | C20orf98    |
| Observed MW: | 26 kDa      |
| UniprotID :  | Q9GZP1      |

### Application


|                       |                                   |
|-----------------------|-----------------------------------|
| Reactivity:           | Human                             |
| Tested Application:   | ELISA, WB, IHC                    |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image:                |                                   |

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human brain using FNab05860(NRSN2 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab05860(NRSN2 antibody) at dilution of 1:500