

## anti- NUP205 antibody

### Product Information

Catalog No.:	FNab05923
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Plays a role in the nuclear pore complex(NPC) assembly and/or maintenance. May anchor NUP62 and other nucleoporins, but not NUP153 and TPR, to the NPC.

### Immunogen information

Immunogen:	nucleoporin 205kDa
Synonyms:	C7orf14, KIAA0225
Observed MW:	205 kDa
UniprotID :	Q92621

### Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human kidney slide using FNab05923(NUP205 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab05923(NUP205 Antibody) at dilution of 1:2000