

anti- OCIAD1 antibody

Product Information

Catalog No.:	FNab05958
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	OCIA domain containing 1
Synonyms:	OCIA
Observed MW:	29-35 kDa
UniprotID :	Q9NX40

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver using FNab05958(OCIAD1 antibody) at dilution of 1:100



IP Result of anti-OCIAD1 (IP:FNab05958, 4ug; Detection:FNab05958 1:300) with HeLa cells lysate 3200ug.



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab05958(OCIAD1 antibody) at dilution of 1:300