

anti- P2RX5 antibody

Product Information

Catalog No.:	FNab06070
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Receptor for ATP that acts as a ligand-gated ion channel.

Immunogen information

Immunogen:	purinergic receptor P2X, ligand-gated ion channel, 5
Synonyms:	P2X5
Observed MW:	44 kDa,49 kDa,116 kDa
UniprotID :	Q93086

Application


Reactivity:	Human
Tested Application:	ELISA, WB, FC
Recommended dilution:WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab06070(P2RX5 antibody) at dilution of 1:600