

## anti- PAX3 antibody

### Product Information

Catalog No.:	FNab06169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Transcription factor that may regulate cell proliferation, migration and apoptosis. Involved in neural development and myogenesis.

### Immunogen information

Immunogen:	paired box 3
Synonyms:	CDHS, HUP2, paired box 3, Paired box protein Pax 3, PAX3, WS1
Observed MW:	60 kDa
Uniprot ID :	P23760

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-PAX3 (IP: FNab06169, 4ug;  
Detection: FNab06169 1:500) with mouse skin  
tissue lysate 3000ug.



HEK-293 cells were subjected to SDS PAGE  
followed by western blot with FNab06169(PAX3  
antibody) at dilution of 1:1000