

anti- PAX3 antibody

Product Information

Catalog No.:	FNab06170
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Transcription factor that may regulate cell proliferation, migration and apoptosis. Involved in neural development and myogenesis.

Immunogen information

Immunogen:	paired box 3
Synonyms:	CDHS, HUP2, paired box 3, Paired box protein Pax 3, PAX3, WS1
Observed MW:	60 kDa
UniprotID :	P23760

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human lung tissue using FNab06170(PAX3 antibody) at dilution of 1:50



mouse skin tissue were subjected to SDS PAGE followed by western blot with FNab06170(PAX3 antibody) at dilution of 1:500