

anti- PDSS2 antibody

Product Information

Catalog No.:	FNab06292
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Supplies decaprenyl diphosphate, the precursor for the side chain of the isoprenoid quinones ubiquinone-10.

Immunogen information

Immunogen:	prenyl(decaprenyl) diphosphate synthase, subunit 2
Synonyms:	C6orf210, DLP1
Observed MW:	27-30KD,44-47KD
UniprotID :	Q86YH6

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human kidney using FNab06292(PDSS2 antibody) at dilution of 1:50



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab06292(PDSS2 antibody) at dilution of 1:300