

anti- Pepsinogen I antibody

Product Information

Catalog No.:	FNab06308
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

This antibody can recognize PGA3, PGA4, PGA5.

Immunogen information

Immunogen:	pepsinogen 5, group I(pepsinogen A)
Synonyms:	Pepsin A, pepsinogen 5, group I (pepsinogen A), Pepsinogen I, PGA3, PGA4, PGA5
Observed MW:	42 kDa
UniprotID :	P0DJJ9

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human stomach cancer tissue slide using FNab06308(PGA5 Antibody) at dilution of 1:100



SGC-7901 cells were subjected to SDS PAGE followed by western blot with FNab06308(PGA5 antibody) at dilution of 1:1000