

anti- PEX5L antibody

Product Information

Catalog No.:	FNab06332
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Accessory subunit of hyperpolarization-activated cyclic nucleotide-gated(HCN) channels, regulating their cell-surface expression and cyclic nucleotide dependence.

Immunogen information

Immunogen:	peroxisomal biogenesis factor 5-like
Synonyms:	Peroxin 5 related protein, PEX2 related protein, PEX5 like protein, PEX5 related protein, PEX5L, PEX5R, Pex5Rp, PXR2, PXR2B
Observed MW:	68 kDa
UniprotID :	Q8IYB4

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab06332(PEX5L Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab06332(PEX5L antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com