

anti- PFKFB3-Specific antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab06339 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|--|
| Immunogen: | 6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 3 |
| Synonyms: | 6PF 2 K/Fru 2,6 P2ase 3, FLJ37326, IPFK2, PFK/FBPase 3, PFK2, PFKFB3 |
| Observed MW: | 58kDa |
| Uniprot ID : | Q16875 |

Application


| | |
|-----------------------|---|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, IHC, WB, IF |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:50-1:200; IF: 1:20-1:200 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab06339(PFKFB3 Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab06339(PFKFB3-Specific antibody) at dilution of 1:600

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com