

## anti- PHACTR1 antibody

### Product Information

Catalog No.:	FNab06370
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Binds actin monomers(G actin) and plays a role in the reorganization of the actin cytoskeleton and in formation of actin stress fibers. Plays a role in cell motility.Plays a role in the formation of tubules by endothelial cells.Regulates PPP1CA activity.Required for normal cell survival.

### Immunogen information

Immunogen:	phosphatase and actin regulator 1
Synonyms:	KIAA1733, RPEL1
Observed MW:	66-70 kDa
UniprotID :	Q9C0D0

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab06370( PHACTR1 Antibody) at dilution of 1:2000

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)