

## anti- PIGP antibody

### Product Information

Catalog No.:	FNab06444
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Part of the complex catalyzing the transfer of N-acetylglucosamine from UDP-N-acetylglucosamine to phosphatidylinositol, the first step of GPI biosynthesis.

### Immunogen information

Immunogen:	phosphatidylinositol glycan anchor biosynthesis, class P
Synonyms:	DCRC, DSCR5, DSCRC
Observed MW:	None
Uniprot ID :	P57054

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA
Recommended dilution:	None
Image:	

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)