

## anti- PJA1 antibody

### Product Information

Catalog No.:	FNab06477
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Has E2-dependent E3 ubiquitin-protein ligase activity. Ubiquitinates MAGED1 antigen leading to its subsequent degradation by proteasome(By similarity). May be involved in protein sorting.

### Immunogen information

Immunogen:	praja ring finger 1
Synonyms:	RNF70
Observed MW:	71-80 kDa
UniprotID :	Q8NG27

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-PJA1 (IP:FNab06477, 3ug;  
Detection:FNab06477 1:300) with mouse brain  
tissue lysate 3600ug.



mouse brain tissue were subjected to SDS PAGE  
followed by western blot with FNab06477(PJA1  
antibody) at dilution of 1:200