

## anti- PLAA antibody

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNab06508   |
| Size:         | 100µg   |
| Form:         | liquid  |
| Purification: | Immunogen affinity purified   |
| Purity:       | ≥95% as determined by SDS-PAGE  |
| Host:         | Rabbit  |
| Clonality:    | polyclonal  |
| Clone ID:     | None  |
| IsoType:      | IgG   |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

### Background

Involved in the maintenance of ubiquitin levels.

### Immunogen information

|              |                                     |
|--------------|-------------------------------------|
| Immunogen:   | phospholipase A2-activating protein |
| Synonyms:    | PLAP                                |
| Observed MW: | 72-74 kDa                           |
| UniprotID :  | Q9Y263                              |

### Application


|                       |                                    |
|-----------------------|------------------------------------|
| Reactivity:           | Human                              |
| Tested Application:   | ELISA, WB, IP                      |
| Recommended dilution: | WB: 1:500-1:2000; IP: 1:200-1:2000 |
| Image:                |                                    |

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



IP Result of anti-PLAA (IP:FNab06508, 4ug;  
Detection:FNab06508 1:600) with K-562 cells  
lysate 3200ug.



K-562 cells were subjected to SDS PAGE  
followed by western blot with FNab06508(PLAA  
antibody) at dilution of 1:500