

anti- phospho-PLK1(S326) antibody

Product Information

Catalog No.:	FNab06547
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	polo-like kinase 1(Drosophila)
Synonyms:	None
Observed MW:	60 kDa
Uniprot ID :	P53350

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon using FNab06547(PLK1-phospho antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab06547(PLK1-phospho(S326) antibody) at dilution of 1:200