

## anti- PLP2 antibody

### Product Information

Catalog No.:	FNab06555
Size:	100μ g
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May play a role in cell differentiation in the intestinal epithelium.

### Immunogen information

Immunogen:	proteolipid protein 2
Synonyms:	A4, A4 LSB, Intestinal membrane A4 protein, PLP2, Proteolipid protein 2
Calculated MW:	15-17kDa
Uniprot ID :	Q04941

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human colon using FNab06555(PLP2 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06555(PLP2 antibody) at dilution of 1:500