

anti- POLR2B-Specific antibody

Product Information

Catalog No.:	FNab06620
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	polymerase(RNA) II(DNA directed) polypeptide B, 140kDa
Synonyms:	None
Observed MW:	None
UniprotID :	P30876

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:5000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain using FNab06620(POLR2B-Specific antibody) at dilution of 1:400



IP Result of anti-POLR2B-Specific (IP:FNab06620, 4ug; Detection:FNab06620 1:300) with HeLa cells lysate 3000ug.

A431 cells were subjected to SDS PAGE followed by western blot with FNab06620(POLR2B-Specific antibody) at dilution of 1:500

