

## anti- PPIL2 antibody

### Product Information

Catalog No.:	FNab06680
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

May catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides thereby assisting the folding of proteins(By similarity). May also function as a chaperone, playing a role in transport to the cell membrane of BSG for instance(PubMed:15946952). May also have a protein ubiquitin ligase activity acting as an E3 ubiquitin protein ligase or as an ubiquitin-ubiquitin ligase promoting elongation of ubiquitin chains on substrates. By mediating 'Lys-48'-linked polyubiquitination of proteins could target them for proteasomaldegradation(PubMed:11435423).

### Immunogen information

Immunogen:	peptidylprolylisomerase(cyclophilin)-like 2
Synonyms:	None
Observed MW:	59 kDa
UniprotID :	Q13356

### Application

Reactivity:	Human, Mouse, Rat
-------------	-------------------

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Tested Application: ELISA, WB

Recommended dilution: WB: 1:500-1:2000

Image:



human kidney tissue were subjected to SDS PAGE followed by western blot with FNab06680(PPIL2 antibody) at dilution of 1:400

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)