

anti- PPP1R2P9 antibody

Product Information

Catalog No.:	FNab06706
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Functions as a protein phosphatase inhibitor. It inhibits activity of the catalytic subunit of PP1 and weakly inhibits the activity of myosin-associated phosphates.

Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 2 pseudogene 9
Synonyms:	None
Observed MW:	35 kDa
Uniprot ID :	O14990

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis using FNab06706(PPP1R2P9 antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab06706(PPP1R2P9 antibody) at dilution of 1:300