

anti- PPP1R9A antibody

Product Information

Catalog No.:	FNab06710
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Binds to actin filaments(F-actin) and shows cross-linking activity. Binds along the sides of the F-actin. May be involved in neurite formation. Inhibits protein phosphatase 1-alpha activity(By similarity).

Immunogen information

Immunogen:	protein phosphatase 1, regulatory(inhibitor) subunit 9A
Synonyms:	KIAA1222
Observed MW:	180 kDa
UniprotID :	Q9ULJ8

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab06710(PPP1R9A antibody) at dilution of 1:300