

anti- PPP2R3A antibody

Product Information

Catalog No.:	FNab06719
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

The B regulatory subunit might modulate substrate selectivity and catalytic activity, and also might direct the localization of the catalytic enzyme to a particular subcellular compartment.

Immunogen information

Immunogen:	protein phosphatase 2(formerly 2A), regulatory subunit B", alpha
Synonyms:	PPP2R3
Observed MW:	130 kDa
UniprotID :	Q06190

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human brain using FNab06719(PPP2R3A antibody) at dilution of 1:100



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab06719(PPP2R3A antibody) at dilution of 1:300