

## anti- PPRC1 antibody

### Product Information

Catalog No.:	FNab06738
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Acts as a coactivator during transcriptional activation of nuclear genes related to mitochondrial biogenesis and cell growth. Involved in the transcription coactivation of CREB and NRF1 target genes.

### Immunogen information

Immunogen:	peroxisome proliferator-activated receptor gamma, coactivator-related 1
Synonyms:	KIAA0595, PGC 1 related coactivator, PPRC1, PRC, RP11 302K17.6
Observed MW:	177 kDa
Uniprot ID :	Q5VV67

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200 - 1:2000
Image:	

1


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab06738( PPRC1 Antibody) at dilution of 1:300