

anti- PRPH2 antibody

Product Information

Catalog No.:	FNab06822
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

May function as an adhesion molecule involved in stabilization and compaction of outer segment disks or in the maintenance of the curvature of the rim. It is essential for disk morphogenesis.

Immunogen information

Immunogen:	peripherin 2(retinal degeneration, slow)
Synonyms:	AOFMD, AVMD, Peripherin 2, Peripherin2, PRPH, PRPH2, rd2, RDS, RP7, Tetraspanin 22, Tspan 22, TSPAN22
Observed MW:	39 kDa
Uniprot ID :	P23942

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab06822(PRPH2 antibody) at dilution of 1:1000