

anti- PRSS2 antibody

Product Information

Catalog No.:	FNab06833
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

In the ileum, may be involved in defensin processing, including DEFA5.

Immunogen information

Immunogen:	protease, serine, 2(trypsin 2)
Synonyms:	Anionic trypsinogen, PRSS2, Serine protease 2, TRY2, TRY8, TRYP2, Trypsin 2, Trypsin II
Observed MW:	26kd,35kd
Uniprot ID :	P07478

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB,IHC
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:200-1:600
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab06833 (PRSS2 Antibody) at dilution of 1:400



Raji cells were subjected to SDS PAGE followed by western blot with FNab06833 (PRSS2 antibody) at dilution of 1:500