

anti- RAB27A antibody

Product Information

Catalog No.: FNab07010

Size: 100μg Form: liquid

Purification: Protein A+G purification

Purity: ≥95% as determined by SDS-PAGE

Host: Mouse

Clonality: monoclonal

Clone ID: 9F7
IsoType: IgG1

Storage: PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12

months (Avoid repeated freeze / thaw cycles.)

Background

Plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

Immunogen information

Immunogen: RAB27A, member RAS oncogene family

Synonyms: GS2, GTP binding protein Ram, HsT18676, Rab 27, RAB27, RAB27A,

RAM, Ras related protein Rab 27A

Observed MW: 25kd Uniprot ID: P51159

Application

Reactivity: Human, Mouse, Rat

Tested Application: ELISA, WB, IHC, IF, IP

Recommended dilution: WB: 1:500-1:2000; IP:1:500-1:2000; IHC: 1:50-1:500; IF: 1:10-1:100

Image:


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)


Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>

1






Immunohistochemistry of paraffin-embedded human kidney using FNab07010(RAB27A antibody) at dilution of 1:50



Immunofluorescent analysis of HepG2 cells, using RAB27A antibody FNab07010 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red)



K-562 cells were subjected to SDS PAGE followed by western blot with FNab07010(RAB27A antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel :(0086)027-87384275 Fax: (0086)027-87800889 <u>www.fn-test.com</u>