

anti- RAB27A-specific antibody

Product Information

Catalog No.:	FNab07011
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	RAB27A, member RAS oncogene family
Synonyms:	GS2, GTP binding protein Ram, HsT18676, Rab 27, RAB27, RAB27A, RAM, Ras related protein Rab 27A
Observed MW:	25kd
Uniprot ID :	P51159

Application

Reactivity:	Human; other species are not tested. Please decide the specificity by homology
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using FNab07011(RAB27A-specific antibody) at dilution of 1:50



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab07011(RAB27A-specific Antibody) at dilution of 1:600