

## anti- RASGEF1A antibody

### Product Information

Catalog No.:	FNab07129
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Guanine nucleotide exchange factor(GEF) with specificity for RAP2A, KRAS, HRAS, and NRAS(in vitro). Plays a role in cell migration.

### Immunogen information

Immunogen:	RasGEF domain family, member 1A
Synonyms:	None
Observed MW:	55 kDa
UniprotID :	Q8N9B8

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human brain using FNab07129(RASGEF1A antibody) at dilution of 1:100



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07129(RASGEF1A Antibody) at dilution of 1:1000