

## anti- RBAP48 antibody

### Product Information

Catalog No.:	FNab07148
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

None

### Immunogen information

Immunogen:	retinoblastoma binding protein 4
Synonyms:	Chromatin assembly factor 1 subunit C,CAF-I 48 kDa subunit,Nucleosome-remodeling factor subunit RBAP48,RBBP-4,Retinoblastoma-binding protein p48
Observed MW:	48-50 kDa
Uniprot ID :	Q09028

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC,IF
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200; IF: 1:50-1:500
Image:	

1


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human cervical cancer using FNab07148(RBBP4 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07148(RBBP4 antibody) at dilution of 1:500