

anti- ROBO1-Specific antibody

Product Information

Catalog No.:	FNab07372
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	roundabout, axon guidance receptor, homolog 1
Synonyms:	Deleted in U twenty twenty, DUTT1, H Robo 1, ROBO1, ROBO1-Specific, Roundabout homolog 1, SAX3
Observed MW:	250 kDa
Uniprot ID :	Q9Y6N7

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07372(ROBO1-Specific antibody) at dilution of 1:1000