

# anti- RPL23 antibody

## **Product Information**

Catalog No.:	FNab07422
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	$\geq$ 95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)
Rackground	

### Background

None

# **Immunogen information**

Immunogen:	ribosomal protein L23
Synonyms:	None
Observed MW:	15 kDa
UniprotID :	P62829

## Application


Reactivity:	Human, Mouse, Rat		
Tested Application:	ELISA, WB, IHC		
Recommended dilution: WB: 1:500-1:2000; IHC: 1:20-1:200			
Image:			

#### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com





Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab07422(RPL23 Antibody) at dilution of 1:200

HeLa cells were subjected to SDS PAGE followed by western blot with FNab07422(RPL23 antibody) at dilution of 1:200



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275 Fax: (0086)027-87800889 www.fn-test.com