

anti- SERPINA10 antibody

Product Information

Catalog No.:	FNab07740
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4F10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Inhibits activity of the coagulation protease factor Xa in the presence of PROZ, calcium and phospholipids. Also inhibits factor XIa in the absence of cofactors.

Immunogen information

Immunogen:	serpin peptidase inhibitor, clade A(alpha-1 antiproteinase, antitrypsin), member 10
Synonyms:	ZPI
Observed MW:	51 kDa
UniprotID :	Q9UK55

Application


Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:1500-1:6000; IF: 1:50-1:500; IHC: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human liver using FNab07740(SERPINA10 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07740(SERPINA10 antibody) at dilution of 1:300