

anti- SGF29 antibody

Product Information

Catalog No.:	FNab07805
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in transcriptional regulation, through association with histone acetyltransferase(HAT) SAGA-type complexes like the TFTC-HAT, ATAC or STAGA complexes. Specifically recognizes and binds methylated 'Lys-4' of histone H3(H3K4me), with a preference for trimethylated form(H3K4me3). In the SAGA-type complexes, required to recruit complexes to H3K4me. May be involved in MYC-mediated oncogenic transformation.

Immunogen information

Immunogen:	coiled-coil domain containing 101
Synonyms:	CCDC101
Observed MW:	33 kDa
UniprotID :	Q96ES7

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab07805 (CCDC101 Antibody) at dilution of 1:400



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07805 (CCDC101 antibody) at dilution of 1:1000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com