

anti- SLC12A2 antibody

Product Information

Catalog No.:	FNab07906
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Electrically silent transporter system. Mediates sodium and chloride reabsorption. Plays a vital role in the regulation of ionic balance and cell volume.

Immunogen information

Immunogen:	solute carrier family 12(sodium/potassium/chloride transporters), member 2
Synonyms:	Basolateral Na K Cl symporter, BSC, BSC2, NKCC1, SLC12A2
Observed MW:	140 kDa
Uniprot ID :	P55011

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IP: 1:200-1:2000; IHC: 1:20-1:200
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human colon cancer using FNab07906(SLC12A2 antibody) at dilution of 1:100



Various lysates were subjected to SDS PAGE followed by western blot with FNab07906(SLC12A2 antibody) at dilution of 1:1000