

anti- SLC44A1 antibody

Product Information

Catalog No.:	FNab07964
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Choline transporter. May be involved in membrane synthesis and myelin production.

Immunogen information

Immunogen:	solute carrier family 44, member 1
Synonyms:	CD92, CDW92, CTL1
Observed MW:	48-55 kDa, 70-73 kDa
UniprotID :	Q8WWI5

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China (430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



U-937 cells were subjected to SDS PAGE followed by western blot with FNab07964(SLC44A1 Antibody) at dilution of 1:1000