

## anti- SOX10 antibody

### Product Information

Catalog No.:	FNab08120
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Transcription factor that seems to function synergistically with the POU domain protein TST-1/OCT6/SCIP. Could confer cell specificity to the function of other transcription factors in developing and mature glia(By similarity).

### Immunogen information

Immunogen:	SRY(sex determining region Y)-box 10
Synonyms:	DOM, SOX10, Transcription factor SOX 10, WS2E, WS4
Observed MW:	50-60 kDa
UniprotID :	P56693

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:500
Image:	

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human brain tissue slide using FNab08120(SOX10 antibody) at dilution of 1:100



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab08120(SOX10 Antibody) at dilution of 1:400