

## anti- SR140 antibody

### Product Information

Catalog No.:	FNab08211
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

#### Immunogen information

Immunogen:	U2-associated SR140 protein
Synonyms:	KIAA0332, SR140, U2 associated protein SR140, U2 associated SR140 protein, U2SURP
Observed MW:	118 kDa
UniprotID :	O15042

#### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:20 - 1:200
Image:	

#### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human breast cancer using FNab08211(SR140 antibody) at dilution of 1:100



Hela cells were subjected to SDS PAGE followed by western blot with FNab08211(SR140 antibody) at dilution of 1:1000