

anti- SRC 1 antibody

Product Information

Catalog No.:	FNab08217
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	nuclear receptor coactivator 1
Synonyms:	bHLHe42, bHLHe74, F SRC 1, KAT13A, NCoA 1, NCOA1, nuclear receptor coactivator 1, Protein Hin 2, RIP160, SRC 1, SRC1, Steroid receptor coactivator 1
Observed MW:	160 kDa
Uniprot ID :	Q15788

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	

1


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human testis using FNab08217(NCOA1 antibody) at dilution of 1:100



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08217(NCOA1 antibody) at dilution of 1:400