

anti- SRD5A1 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab08219 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Converts testosterone into 5-alpha-dihydrotestosterone and progesterone or corticosterone into their corresponding 5-alpha-3-oxosteroids. It plays a central role in sexual differentiation and androgen physiology.

Immunogen information

| | |
|--------------|---|
| Immunogen: | steroid-5-alpha-reductase, alpha polypeptide 1(3-oxo-5 alpha-steroid delta 4-dehydrogenase alpha 1) |
| Synonyms: | None |
| Observed MW: | 29 kDa |
| UniprotID : | P18405 |

Application

| | |
|--------------------------|--------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:500-1:2000 |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

Image:



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08219(SRD5A1 Antibody) at dilution of 1:1000