

anti- ST7L antibody

Product Information

Catalog No.:	FNab08274
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	suppression of tumorigenicity 7 like
Synonyms:	ST7R
Observed MW:	60-66 kDa
UniprotID :	Q8TDW4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

IP Result of anti-ST7L (IP:FNab08274, 3ug;
Detection:FNab08274 1:1500) with mouse testis
tissue lysate 5000ug.



mouse testis tissue were subjected to SDS PAGE
followed by western blot with FNab08274(ST7L
antibody) at dilution of 1:1500