

## anti- STAR antibody

### Product Information

Catalog No.:	FNab08288
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

Plays a key role in steroid hormone synthesis by enhancing the metabolism of cholesterol into pregnenolone. Mediates the transfer of cholesterol from the outer mitochondrial membrane to the inner mitochondrial membrane where it is cleaved to pregnenolone.

### Immunogen information

Immunogen:	steroidogenic acute regulatory protein
Synonyms:	STARD1
Observed MW:	33 kDa
Uniprot ID :	P49675

### Application

Reactivity:	Human, Mouse, Rat, Pig
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	


### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

[www.fn-test.com](http://www.fn-test.com)



Immunohistochemistry of paraffin-embedded human ovary tissue slide using FNab08288(STAR Antibody) at dilution of 1:50



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with FNab08288(STAR antibody) at dilution of 1:500