

anti- STK38L antibody

Product Information

Catalog No.:	FNab08338
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Involved in the regulation of structural processes in differentiating and mature neuronal cells.

Immunogen information

Immunogen:	serine/threonine kinase 38 like
Synonyms:	KIAA0965, NDR2, NDR2 protein kinase, Nuclear Dbf2 related kinase 2, STK38L
Observed MW:	54 kDa
UniprotID :	Q9Y2H1

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded rat heart using FNab08338(STK38L antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08338(STK38L antibody) at dilution of 1:800