

anti- STMN3 antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab08342 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Protein A+G purification |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 6D5 |
| IsoType: | IgG2a |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Exhibits microtubule-destabilizing activity, which is antagonized by STAT3.

Immunogen information

| | |
|--------------|-----------------|
| Immunogen: | stathmin-like 3 |
| Synonyms: | SCLIP |
| Observed MW: | 21 kDa |
| Uniprot ID : | Q9NZ72 |

Application

| | |
|--------------------------|--------------|
| Reactivity: | Human |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:500-1:2000 |
| Image: | |


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



human brain tissue were subjected to SDS PAGE followed by western blot with FNab08342(STMN3 antibody) at dilution of 1:1000