

anti- TAB2 antibody

Product Information

Catalog No.:	FNab08465
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

Adapter linking MAP3K7/TAK1 and TRAF6.Promotes MAP3K7 activation in the IL1 signaling pathway. The binding of 'Lys-63'-linked polyubiquitin chains to TAB2 promotes autophosphorylation of MAP3K7 at 'Thr-187'. Involved in heart development.

Immunogen information

Immunogen:	mitogen-activated protein kinase kinasekinase 7 interacting protein 2
Synonyms:	KIAA0733, MAP3K7IP2
Observed MW:	77 kDa
UniprotID :	Q9NYJ8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IF, WB, FC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	


1

Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275


Fax: (0086)027-87800889 www.fn-test.com



Immunohistochemistry of paraffin-embedded human malignant melanoma using FNab08465(MAP3K7IP2 antibody) at dilution of 1:100



IP Result of anti-TAB2 (IP:FNab08465, 4ug; Detection:FNab08465 1:600) with HepG2 cells lysate 3600ug.



L02 cells were subjected to SDS PAGE followed by western blot with FNab08465(TAB2 antibody) at dilution of 1:2000