

## anti- TAX1BP1 antibody

### Product Information

Catalog No.:	FNab08506
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

TAX1BP1 (also named as T6BP and TXBP151) inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. It is degraded by caspase-3-like family proteins upon TNF-induced apoptosis. TAX1BP1 may also play a role in the pro-inflammatory cytokine IL-1 signaling cascade. TAX1BP1 binds TRAF6 but not TRAFs 1-5. TAX1BP1 has 4 isoforms with MW 65-68 kDa and 86-91 kDa.

### Immunogen information

Immunogen:	Tax1(human T-cell leukemia virus type I) binding protein 1
Synonyms:	CALCOCO3, T6BP, Tax1 binding protein 1, TAX1BP1, TRAF6 binding protein, TXBP151
Observed MW:	68 kDa, 85 kDa
Uniprot ID :	Q86VP1

### Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC

### Wuhan Fine Biotech Co., Ltd.


B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Recommended dilution: WB: 1:500 - 1:2000; IHC: 1:50 - 1:200

Image:



Immunohistochemistry of paraffin-embedded rat kidney tissue slide using FNab08506( TAX1BP1 Antibody) at dilution of 1:100



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab08506(TAX1BP1 antibody) at dilution of 1:500