

anti- TDP-43 antibody

Product Information

Catalog No.:	FNab08576
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4G1
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	TAR DNA binding protein
Synonyms:	ALS10, TAR DNA binding protein, TAR DNA binding protein 43, TARDBP, TDP 43, TDP43
Observed MW:	43 kDa
Uniprot ID :	Q13148

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:5000-1:20000; IHC: 1:1000-1:2000
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)


Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab08576(TDP-43 antibody) at dilution of 1:1000



Hela cells(shcontrol and shRNA of TDP43) were subjected to SDS PAGE followed by western blot with FNab08576(TDP-43 antibody) at dilution of 1:10000