

anti- TNFRSF6B antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab08826 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Background

Decoy receptor that can neutralize the cytotoxic ligands TNFSF14/LIGHT, TNFSF15 and TNFSF6/FASL. Protects against apoptosis.

Immunogen information

| | |
|--------------|--|
| Immunogen: | tumor necrosis factor receptor superfamily, member 6b, decoy |
| Synonyms: | DCR3, TR6 |
| Observed MW: | None |
| UniprotID : | O95407 |

Application


| | |
|--------------------------|--------------|
| Reactivity: | Human |
| Tested Application: | ELISA, IP |
| Recommended dilution:IP: | 1:200-1:2000 |
| Image: | |

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com



IP Result of anti-TNFRSF6B (IP:FNab08826, 4ug;
Detection:FNab08826 1:300) with HeLa cells
lysate 2480ug.