

anti- TPI1 antibody

Product Information

Catalog No.:	FNab08887
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2C10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	triosephosphate isomerase 1
Synonyms:	TPI
Observed MW:	27 kDa
UniprotID :	P60174

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000-1:4000
Image:	

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :(0086)027-87384275

Fax: (0086)027-87800889 www.fn-test.com

HepG2 cells were subjected to SDS PAGE followed by western blot with FNab08887(TPI1 antibody) at dilution of 1:2000

