

## anti- Transgelin-2-specific antibody

### Product Information

Catalog No.:	FNab08937
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background

The transgelin family is a group of proteins that belong to 22kd actin-related corpnin superfamily. Of all three isoforms, transgelin 1 is the best characterized. Transgelin 1, also known as SM22alpha, is a specific marker for differentiated smooth muscle cells. Transgelin 2, also known as SM22 beta, is expressed by both smooth muscle and non-smooth muscle cells in a temporally and spatially regulated pattern. Trangenlin 3, also known as NP25, is only found in highly differentiated neuronal cells. This antibody can bind transgelin 2 specifically, but not trangelin 1 or 3.

### Immunogen information

Immunogen:	transgelin 2
Synonyms:	Epididymis tissue protein Li 7e,SM22-alpha homolog,TAGLN2,KIAA0120,CDABP0035
Observed MW:	22 kDa
UniprotID :	P37802

### Application

### Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve.East Lake High-Tech Development Zone.Wuhan, Hubei, China(430206)

Tel :( 0086)027-87384275


Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

Reactivity: Human, Mouse, Rat


Tested Application: ELISA, IHC, IF, IP, FC, WB

Recommended dilution: WB: 1:200-1:2000; IP: 1:200-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100; FC:N/A

Image:



Immunohistochemistry of paraffin-embedded human kidney using FNab08937(Transgelin-2-specific antibody) at dilution of 1:50



IP Result of anti-Transgelin-2-specific (IP:FNab08937, 3ug; Detection:FNab08937 1:500) with L02 cells lysate 1000ug.

**Wuhan Fine Biotech Co., Ltd.**

B9 Bld, High-Tech Medical Devices Park, No. 818 GaoxinAve. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889 [www.fn-test.com](http://www.fn-test.com)

rat stomach tissue were subjected to SDS PAGE followed by western blot with FNab08937(Transgelin-2-specific Antibody) at dilution of 1:600

